

KEY FEATURES

- Smooth, 3-Step or 5-Step Turret Magnification System
- LED Illumination System
- APOchromatic Lenses
- 175-400mm Working Distance Objective Lens
- Digital, Live Video and/or HD Video Capture Card System
- Widefield Eyepieces
- Straight Binocular Head (Optional: Fixed Inclined or Inclinable)
- Green Filter
- Lifetime Warranty on Optics and Mechanics

COLPOSCOPES

FOCUS ON PRECISION

SEILER 900 SERIES LED OSCOPES





935 3-Step Colposcope
Small Diameter Base on Roller Casters

with Adjustable Vertical Post Height 3 Levels of Magnification



955 5-Step Colposcope

Comfort Mounting on Swing Arm Base Provides Small Footprint 5 Levels of Magnification



985 OTS Colposcope

Over-the-Shoulder Mounting Removes Obstructions from the Working Field 5 Levels of Magnification

Live Capture of HD Images and Video Available on All Models



0-550° BINOCULAR HEAD
The inclinable head allows for a full 0-220° rotation. Seiler offers: 0-220° inclinable,



DIGITAL DSLR CAMERA



<u>ASSISTANT</u> **TEACHING HEAD**



BEAMSPLITTER

Specially designed to create the highest quality imagery. The port contains a 50/50 side for digital and 80/20 for ideo. Optional 50/50-50/50



VARIO FOCUS Never change an objective again. Stay in focus from 200mm to 350mm by simply turning the fine focus knob.



GO PRO HERO 5



FIXED 45° HEAD FOR HRA PROCEDURES

NEW HD VIDEO CAPTURE CARD SYSTEM





