

ACCU-CAM 300 Eyepiece Camera

The ACCU-CAM 300 Eyepiece Camera is the perfect solution for those on a tight budget who need to capture images and video!



CAT# AU-300-EP

- 3.0 MP CMOS eyepiece camera
- Sensor size: 1/2.7"
- USB 2.0 direct connect to your PC!
- Includes our new, feature-rich AccuView software for image analysis, annotation, measurement, white balance & more!
- Captures digital still and video
- Windows compatible
- Easy to use
- Excellent color reproduction
- Fits into a 23.2mm I.D. eye tube or photo tube
- Optional adapter (CAT #5036) available for 30mm I.D. eyepieces
- What's in the box: ACCU-CAM 300 eyepiece camera, USB 2.0 mini cable, AccuView software CD



*Delivering quality
through vision*



73 Mall Drive • Commack, NY 11725
631-864-1000 (P) • 631-543-8900 (F)
email: info@accu-scope.com • www.accu-scope.com